## SARS-CoV-2-Antigenschnelltest

### Packungsbeilage

REF L031-11815 Deutsch

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BIBLIOGRAPHIE

- Shuo Su, Gary Wong, Weifeng Shi, et al. Epidemiology, Genetic recombination, and pathogenesis of coronaviruses. Trends in Microbiology, Juni 2016, Band 24, Nr. 6: 490-502
- Susan R. Weiss, Julian L. Leibowitz, Coronavirus Pathogenesis, Advances in Virus Research, Band 81: 85-164

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**ACON Biotech (Hangzhou) Co., Ltd.**  
No.210 Zhenzhong Road, West Lake District, Hangzhou, P.R.China, 310030

EC

REP

MedNet EC-REP GmbH  
Borkstrasse 10  
48163 Muenster, Germany

| Method                                                        |          | RT-PCR (Nasopharyngeal Swab Specimens) |          | Total Results |
|---------------------------------------------------------------|----------|----------------------------------------|----------|---------------|
| SARS-CoV-2 Antigen Rapid Test (Nasopharyngeal Swab Specimens) | Results  | Negative                               | Positive |               |
|                                                               | Negative | 175                                    | 3        |               |
|                                                               | Positive | 1                                      | 120      | 121           |
| Total Results                                                 |          | 176                                    | 123      | 299           |

Relative Sensitivity: 97.6% (92.8% - 99.5%)\*      Relative Specificity: 99.4% (96.5% - 99.9%)\*  
Accuracy: 98.7% (96.5% - 99.6%)\*      \*95% Confidence Intervals

Limit of Detection (LOD)

The LOD of SARS-CoV-2 Antigen Rapid Test was established using limiting dilutions of an inactivated viral sample. The viral sample was spiked with negative human nasal and nasopharyngeal sample pool into a series of concentrations. Each level was tested for 30 replicates. The results show that the LOD is 1.6\*10<sup>2</sup> TCID<sub>50</sub>/mL.

Cross-Reactivity (Analytical Specificity) and Microbial Interference

Cross-reactivity was evaluated by testing a panel of related pathogens and microorganisms that are likely to be present in the nasal cavity. Each organism and virus were tested in the absence or presence of heat-inactivated SARS-CoV-2 virus at low positive level.

No cross-reactivity or interference was observed with the following microorganisms when tested at the concentration presented in the table below. The SARS-CoV-2 Antigen Rapid Test does not differentiate between SARS-CoV and SARS-CoV-2.

| Potential Cross-Reactant |                             | Test Concentration                            | Cross-Reactivity<br>(in the absence of SARS-CoV-2 virus) | Interference<br>(in the presence of SARS-CoV-2 virus) |
|--------------------------|-----------------------------|-----------------------------------------------|----------------------------------------------------------|-------------------------------------------------------|
| Virus                    | Adenovirus                  | 1.14 x 10 <sup>6</sup> TCID <sub>50</sub> /mL | No<br>3/3 negative                                       | No<br>3/3 positive                                    |
|                          | Enterovirus                 | 9.50 x 10 <sup>5</sup> TCID <sub>50</sub> /mL | No<br>3/3 negative                                       | No<br>3/3 positive                                    |
|                          | Human coronavirus 229E      | 1.04 x 10 <sup>5</sup> TCID <sub>50</sub> /mL | No<br>3/3 negative                                       | No<br>3/3 positive                                    |
|                          | Human coronavirus OC43      | 2.63 x 10 <sup>5</sup> TCID <sub>50</sub> /mL | No<br>3/3 negative                                       | No<br>3/3 positive                                    |
|                          | Human coronavirus NL63      | 1.0 x 10 <sup>5</sup> TCID <sub>50</sub> /mL  | No<br>3/3 negative                                       | No<br>3/3 positive                                    |
|                          | Human Metapneumovirus       | 1.25 x 10 <sup>5</sup> TCID <sub>50</sub> /mL | No<br>3/3 negative                                       | No<br>3/3 positive                                    |
|                          | MERS-coronavirus            | 7.90 x 10 <sup>5</sup> TCID <sub>50</sub> /mL | No<br>3/3 negative                                       | No<br>3/3 positive                                    |
|                          | Influenza A                 | 1.04 x 10 <sup>5</sup> TCID <sub>50</sub> /mL | No<br>3/3 negative                                       | No<br>3/3 positive                                    |
|                          | Influenza B                 | 1.04 x 10 <sup>5</sup> TCID <sub>50</sub> /mL | No<br>3/3 negative                                       | No<br>3/3 positive                                    |
|                          | Parainfluenza virus 1       | 1.25 x 10 <sup>5</sup> TCID <sub>50</sub> /mL | No<br>3/3 negative                                       | No<br>3/3 positive                                    |
|                          | Parainfluenza virus 2       | 3.78 x 10 <sup>5</sup> TCID <sub>50</sub> /mL | No<br>3/3 negative                                       | No<br>3/3 positive                                    |
|                          | Parainfluenza virus 3       | 1.0 x 10 <sup>5</sup> TCID <sub>50</sub> /mL  | No<br>3/3 negative                                       | No<br>3/3 positive                                    |
|                          | Parainfluenza virus 4       | 2.88 x 10 <sup>6</sup> TCID <sub>50</sub> /mL | No<br>3/3 negative                                       | No<br>3/3 positive                                    |
|                          | Respiratory syncytial virus | 3.15 x 10 <sup>5</sup> TCID <sub>50</sub> /mL | No<br>3/3 negative                                       | No<br>3/3 positive                                    |
|                          | Rhinovirus                  | 3.15 x 10 <sup>5</sup> TCID <sub>50</sub> /mL | No<br>3/3 negative                                       | No<br>3/3 positive                                    |
|                          | Human coronavirus-HKU1      | 1 x 10 <sup>5</sup> copies/mL                 | No<br>3/3 negative                                       | No<br>3/3 positive                                    |
| Bacteria                 | Bordetella pertussis        | 2.83 x 10 <sup>9</sup> CFU/mL                 | No<br>3/3 negative                                       | No<br>3/3 positive                                    |
|                          | Chlamydia trachomatis       | 3.13 x 10 <sup>8</sup> CFU/mL                 | No<br>3/3 negative                                       | No<br>3/3 positive                                    |
|                          | Haemophilus influenza       | 1.36 x 10 <sup>8</sup> CFU/mL                 | No<br>3/3 negative                                       | No<br>3/3 positive                                    |
|                          | Legionella pneumophila      | 4.08 x 10 <sup>9</sup> CFU/mL                 | No<br>3/3 negative                                       | No<br>3/3 positive                                    |
|                          | Mycobacterium tuberculosis  | 1.72 x 10 <sup>7</sup> CFU/mL                 | No<br>3/3 negative                                       | No<br>3/3 positive                                    |
|                          | Mycoplasma pneumoniae       | 7.90 x 10 <sup>7</sup> CFU/mL                 | No<br>3/3 negative                                       | No<br>3/3 positive                                    |

|                         |                                      |                               |                    |                    |
|-------------------------|--------------------------------------|-------------------------------|--------------------|--------------------|
|                         | Staphylococcus aureus                | 1.38 x 10 <sup>7</sup> CFU/mL | No<br>3/3 negative | No<br>3/3 positive |
|                         | Staphylococcus epidermidis           | 2.32 x 10 <sup>9</sup> CFU/mL | No<br>3/3 negative | No<br>3/3 positive |
|                         | Streptococcus pneumoniae             | 1.04 x 10 <sup>8</sup> CFU/mL | No<br>3/3 negative | No<br>3/3 positive |
|                         | Streptococcus pyogenes               | 4.10 x 10 <sup>6</sup> CFU/mL | No<br>3/3 negative | No<br>3/3 positive |
|                         | Pneumocystis jirovecii-S. cerevisiae | 8.63 x 10 <sup>7</sup> CFU/mL | No<br>3/3 negative | No<br>3/3 positive |
|                         | Pseudomonas aeruginosa               | 1.87 x 10 <sup>8</sup> CFU/mL | No<br>3/3 negative | No<br>3/3 positive |
|                         | Chlamydia pneumoniae                 | 1×10 <sup>6</sup> IFU/ml      | No<br>3/3 negative | No<br>3/3 positive |
| Yeast                   | Candida albicans                     | 1.57 x 10 <sup>8</sup> CFU/mL | No<br>3/3 negative | No<br>3/3 positive |
| Pooled human nasal wash |                                      |                               | No<br>3/3 negative | No<br>3/3 positive |

Interfering Substances

The following substances, naturally present in respiratory specimens or that may be artificially introduced into the nasal cavity or nasopharynx, were evaluated. Each substance was tested in the absence or presence of SARS-CoV-2 virus at low positive level. The final concentration of the substances tested are listed below and were found not to affect test performance.

| Interfering Substance                           | Active Ingredient                             | Concentration | Results<br>(in the absence of SARS-CoV-2 virus) | Results<br>(in the presence of SARS-CoV-2 virus) |
|-------------------------------------------------|-----------------------------------------------|---------------|-------------------------------------------------|--------------------------------------------------|
| Endogenous                                      | Biotin                                        | 2.4 mg/mL     | 3/3 negative                                    | 3/3 positive                                     |
|                                                 | Mucin                                         | 0.5% w/v      | 3/3 negative                                    | 3/3 positive                                     |
|                                                 | Whole Blood                                   | 4% v/v        | 3/3 negative                                    | 3/3 positive                                     |
| Afrin Original Nasal Spray                      | Oxymetazoline                                 | 15% v/v       | 3/3 negative                                    | 3/3 positive                                     |
| ALKALOL Allergy Relief Nasal Spray              | Homeopathic                                   | 1:10 Dilution | 3/3 negative                                    | 3/3 positive                                     |
| Chloraseptic Max Sore Throat Lozenges           | Menthol, Benzocaine                           | 1.5 mg/mL     | 3/3 negative                                    | 3/3 positive                                     |
| CVS Health Fluticasone Propionate Nasal Spray   | Fluticasone propionate                        | 5% v/v        | 3/3 negative                                    | 3/3 positive                                     |
| Equate Fast-Acting Nasal Spray                  | Phenylephrine                                 | 15% v/v       | 3/3 negative                                    | 3/3 positive                                     |
| Equate Sore Throat Phenol Oral Anesthetic Spray | Phenol                                        | 15% v/v       | 3/3 negative                                    | 3/3 positive                                     |
| Original Extra Strong Menthol Cough Lozenges    | Menthol                                       | 1.5 mg/mL     | 3/3 negative                                    | 3/3 positive                                     |
| NasalCrom Nasal Spray                           | Cromolyn                                      | 15% v/v       | 3/3 negative                                    | 3/3 positive                                     |
| NeilMed NasoGel for Dry Noses                   | Sodium Hyaluronate                            | 5% v/v        | 3/3 negative                                    | 3/3 positive                                     |
| Throat Lozenge                                  | Dyclonine Hydrochloride                       | 1.5mg/mL      | 3/3 negative                                    | 3/3 positive                                     |
| Zicam Cold Remedy                               | Galphimia glauca, Luffa operculata, Sabadilla | 5% v/v        | 3/3 negative                                    | 3/3 positive                                     |
| Antibiotic                                      | Mupirocin                                     | 10 mg/mL      | 3/3 negative                                    | 3/3 positive                                     |
| Tamiflu                                         | Oseltamivir Phosphate                         | 5 mg/mL       | 3/3 negative                                    | 3/3 positive                                     |
| Antibiotic                                      | Tobramycin                                    | 4 µg/mL       | 3/3 negative                                    | 3/3 positive                                     |
| Mometasone Furoate Nasal Spray                  | Mometasone Furoate                            | 5%v/v         | 3/3 negative                                    | 3/3 positive                                     |

|                                      |      |        |              |              |
|--------------------------------------|------|--------|--------------|--------------|
| Physiological Seawater Nasal Cleaner | NaCl | 15%v/v | 3/3 negative | 3/3 positive |
|--------------------------------------|------|--------|--------------|--------------|

PRECISION

Intra-Assay

Within-run precision was determined using 60 replicates of specimens: negative specimens and SARS-CoV-2 antigen positive specimens. The specimens were correctly identified >99% of the time.

Inter-Assay

Between-run precision was determined using 60 independent assays on the same specimen: negative specimen and SARS-CoV-2 antigen positive specimen. Three different lots of the SARS-CoV-2 Antigen Rapid Test were tested using these specimens. The specimens were correctly identified >99% of the time.

BIBLIOGRAPHY

- Shuo Su, Gary Wong, Weifeng Shi, et al. Epidemiology, Genetic recombination, and pathogenesis of coronaviruses. Trends in Microbiology, June 2016, vol. 24, No. 6: 490-502
- Susan R. Weiss, Julian L. Leibowitz, Coronavirus Pathogenesis, Advances in Virus Research, Volume 81: 85-164

Index of Symbols

|                                                                                     |                                                                                     |                                                                                     |                                   |                                                                                     |                     |
|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|-----------------------------------|-------------------------------------------------------------------------------------|---------------------|
|  | Manufacturer                                                                        |  | Contains sufficient for <n> tests |  | Temperature limit   |
|  | In vitro diagnostic medical device                                                  |  | Use-by date                       |  | Do not reuse        |
|  | Consult instructions for use                                                        |  | Batch code                        |  | Catalogue number    |
|  |  | Authorized representative in the European Community                                 |                                   |  | Date of manufacture |

Index of Contents

|                                      |                               |
|--------------------------------------|-------------------------------|
| <b>SARS-CoV-2 Antigen</b>            | SARS-CoV-2 Antigen            |
| <b>Negative Control Swab</b>         | Negative Control Swab         |
| <b>Positive Control Swab</b>         | Positive Control Swab         |
| <b>Extraction Buffer Tubes</b>       | Extraction Buffer Tubes       |
| <b>Disposable Swabs</b>              | Disposable Swabs              |
| <b>Nasal Swabs</b>                   | Nasal swabs                   |
| <b>Nasopharyngeal Swabs</b>          | Nasopharyngeal Swabs          |
| <b>SARS-CoV-2 Antigen Rapid Test</b> | SARS-CoV-2 Antigen Rapid Test |

**ACON Biotech (Hangzhou) Co., Ltd.**  
No.210 Zhenzhong Road, West Lake District, Hangzhou, P.R.China, 310030

**MedNet EC-REP GmbH**  
Borkstrasse 10  
48163 Muenster, Germany