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BIBLIOGRAPHIE

1. Clémence Magnard, Martine Valette, Michèle Aymard, Bruno Lina. Comparison of Two Nested PCR, Cell Culture, and Antigen Detection for the Diagnosis of Upper Respiratory Tract Infections due to Influenza Viruses. Journal of Medical Virology, 59:215–220 (1999)
2. Shuo Su, Gary Wong, Weifeng Shi, et al. Epidemiology, Genetic recombination, and pathogenesis of coronaviruses. Trends in Microbiology, June 2016, vol. 24, No. 6: 490-502
3. Susan R. Weiss, Julian L. Leibowitz, Coronavirus Pathogenesis, Advances in Virus Research, Volume 81: 85-164

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**ACON Biotech (Hangzhou) Co., Ltd.**
No.210 Zhenzhong Road, West Lake
District, Hangzhou, P.R.China, 310030

ECREP

MedNet EC-REP GmbH
Borkstrasse 10
48163 Muenster, Germany

INTERPRETATION OF RESULTS

(Please refer to the illustration above)

NEGATIVE: Only the colored control lines appear in the control region (C). No colored lines appear in the test line regions (T/A/B). This means that no SARS-CoV-2 and/or Influenza A/B antigen were detected.

COVID 19 POSITIVE:* Two distinct colored lines appear in the SARS-CoV-2 side. One line in the control line region (C) and the other line-in the test line region (T). This means that the presence of SARS-CoV-2 antigen was detected.

Flu A POSITIVE:* Two distinct colored lines appear in the Flu A/B side. One line in the control line region (C) and the other line in the test line region (A). This means that the presence of Influenza A antigen was detected.

Flu B POSITIVE:* Two distinct colored lines appear in the Flu A/B side. One line in the control line region (C) and the other line in the test line region (B). This means that the presence of Influenza B antigen was detected.

Flu A/B POSITIVE:* Three distinct colored lines appear in the Flu A/B side. One line in the control line region (C) and the other two lines in the test line region (A and B). This means that the presence of Influenza A/B antigen was detected.

***NOTE:** The intensity of the color in the test line (T/A/B) may vary depending on the level of antigen present in the specimen. Therefore, any shade of color in the test line regions (T/A/B) should be considered positive.

INVALID RESULT: A control line fails to appear in either the SARS-CoV-2, or the Flu A/B side of the cassette. Insufficient specimen volume or incorrect operation are the most likely reasons for a control line’s failure to appear. Review the procedure and repeat the test with a new test cassette. If the problem persists, discontinue using the test kit immediately and contact your local distributor.

QUALITY CONTROL

Internal procedural controls are included in the test. A colored line appearing in the control line region (C) is an internal procedural control. It confirms sufficient specimen volume was applied and correct procedural technique was followed. Positive and Negative control swabs are supplied with each kit. These control swabs should be used to ensure that the test cassette and that the test procedure is performed correctly. Follow the **“DIRECTIONS FOR USE”** section to perform the control test.

The control swabs can be tested under any of the following circumstances:

- 1. When new lot of tests are used and/or when a new operator performs the test.
- 2. At periodic intervals as dictated by local requirements, and/or by the user’s Quality Control procedures.

Note: The control swab should be inserted directly into the extraction buffer tube. Do not insert the control swab into nostril.

LIMITATIONS

- 1. The SARS-CoV-2 & Influenza A/B Ag Combo Rapid Test should be used for the qualitative detection of SARS-CoV-2, Influenza virus type A and type B nucleoprotein antigens in nasal swab and nasopharyngeal swab only. The intensity of the test line does not necessarily correlate to SARS-CoV-2 or Influenza A/B viral titer in the specimen.
- 2. Specimens should be tested as quickly as possible after specimen collection or at most within one (1) hour after specimen collection, if stored at room temperature (15 - 30°C).
- 3. Use of viral transport media may result in decreased test sensitivity.
- 4. A false-negative test result may happen if the level of antigen in a sample is below the detection limit of the test or if the sample was improperly collected.
- 5. Test results should be correlated with other clinical data available to the physician.
- 6. A positive test result does not rule out co-infections with other pathogens.
- 7. A positive test result of SARS-CoV-2 does not differentiate between SARS-CoV and SARS-CoV-2.
- 8. A negative test result is not intended to rule out other viral or bacterial infections.

PERFORMANCE CHARACTERISTICS

Clinical Sensitivity, Specificity and Accuracy

Nasal Swab Specimens

The clinical performance of SARS-CoV-2 & Influenza A/B Ag Combo Rapid Test for nasal swab specimen was established with 580 nasal swabs collected from individual symptomatic patients who were suspected of COVID-19 or Influenza A/B. The results show that the relative sensitivity and the relative specificity are as follows:

Method		RT-PCR (Nasopharyngeal Swab Specimens)		
		Negative	Positive	Total
SARS-CoV-2 Test Results (Nasal Swab Specimens)	Negative	429	4	433
	Positive	1	146	147
	Total	430	150	580

Relative sensitivity: 97.33% (93.11%-99.19%)* Relative specificity: 99.77% (98.56%-99.99%)*
Accuracy: 99.14% (97.94%-99.69%)* *95% Confidence Intervals

Method		Comparator (Nasopharyngeal Swab Specimens)		
		Negative	Positive	Total
Flu A Test Results (Nasal Swab Specimens)	Negative	503	2	505
	Positive	2	73	75
	Total	505	75	580

Relative sensitivity: 97.33% (90.23%-99.83%)* Relative specificity: 99.60% (98.47%-99.99%)*
Accuracy: 99.31% (98.17%-99.80%)* *95% Confidence Intervals

Method		Comparator (Nasopharyngeal Swab Specimens)		
		Negative	Positive	Total
Flu B Test Results (Nasal Swab Specimens)	Negative	504	3	507
	Positive	1	72	73
	Total	505	75	580

Relative sensitivity: 96.00% (88.42%-99.10%)* Relative specificity: 99.80% (98.77%-99.99%)*
Accuracy: 99.31% (98.17%-99.80%)* *95% Confidence Intervals

Nasopharyngeal Swab Specimens

The clinical performance of SARS-CoV-2 & Influenza A/B Ag Combo Rapid Test for nasopharyngeal swab specimen was established with 395 nasopharyngeal swabs collected from individual symptomatic patients who were suspected of COVID-19 or Influenza A/B. The results show that the relative sensitivity and the relative specificity are as follows:

Method		RT-PCR (Nasopharyngeal Swab Specimens)		
		Negative	Positive	Total
SARS-CoV-2 Test Results (Nasopharyngeal Swab Specimens)	Negative	256	3	259
	Positive	1	135	136
	Total	257	138	395

Relative sensitivity: 97.83% (93.52%-99.54%)* Relative specificity: 99.61% (97.60%-99.99%)*
Accuracy: 98.99% (97.33%-99.70%)* *95% Confidence Intervals

Method		Comparator (Nasopharyngeal Swab Specimens)		
		Negative	Positive	Total
Flu A Test Results (Nasopharyngeal Swab Specimens)	Negative	333	1	334
	Positive	2	59	61
	Total	335	60	395

Relative sensitivity: 98.33% (90.30%-99.99%)* Relative specificity: 99.40% (97.70%-99.98%)*
Accuracy: 99.24% (97.68%-99.85%)* *95% Confidence Intervals

Method		Comparator (Nasopharyngeal Swab Specimens)		
		Negative	Positive	Total
Flu B Test Results (Nasopharyngeal Swab Specimens)	Negative	334	1	335
	Positive	1	59	60
	Total	335	60	395

Relative sensitivity: 98.33% (90.30%-99.99%)* Relative specificity: 99.70% (98.15%-99.99%)*
Accuracy: 99.49% (98.05%-99.99%)* *95% Confidence Intervals

Limit of Detection (LOD)

A SARS-CoV-2 viral sample was spiked with negative human nasal and nasopharyngeal sample pool into a series of concentrations. Each level was tested for 30 replicates. The results show that the LOD of SARS-CoV-2 Antigen test is 1.6*10² TCID₅₀/mL.

The Influenza A and Influenza B virus sample was spiked with negative human nasal and nasopharyngeal sample pool into a series of concentrations. Each level was tested for 30 replicates. The results show that the LOD of the Influenza A/B Antigen test is 6.88*10² TCID₅₀/mL on Influenza A virus and 1.88*10² TCID₅₀/mL on Influenza B virus.

Cross-Reactivity and Interference

No cross-reactivity or interference was observed with the following microorganisms: Adenovirus, Enterovirus, Coronavirus-229E, Coronavirus-NL63, Coronavirus-OC43, Human metapneumovirus, MERS-coronavirus, Parainfluenza virus 1, Parainfluenza virus 2,

Parainfluenza virus 3, Parainfluenza virus 4, Respiratory syncytial virus, Rhinovirus, Human coronavirus- HKU1, Bordetella pertussis, Chlamydia trachomatis, Haemophilus influenza, Legionella pneumophila, Mycobacterium tuberculosis, Mycoplasma pneumonia, Staphylococcus aureus, Staphylococcus epidermidis, Streptococcus pneumonia, Streptococcus pyogenes, Pneumocystis jirovecii-S. cerevisiae, Pseudomonas aeruginosa, Chlamydia pneumonia, Candida albicans.

The interfering substances including Whole Blood, Dafenlin Oxymetazoline Hydrochloride Spray, Mometasone Furoate Nasal Spray, Fluticasone Propionate andPhysiological Seawater Nasal Cleaner with a certain concentration have no interference on the test of SARS-CoV-2 & Influenza A/B Ag Combo Rapid Test.

Precision and Reproducibility

Intra-Assay

Within-run precision was determined using 60 replicates of specimens: negative specimen and SARS-CoV-2/Influenza A/Influenza B antigen positive specimen. The specimens were correctly identified 100% of the time.

Inter-Assay

Between-run precision was determined using 60 independent assays on the same specimen: negative specimen and SARS-CoV-2/influenza A/Influenza B antigen positive specimen. Three different lots of the SARS-CoV-2 & Influenza A/B Ag Combo Rapid Test were tested using these specimens. The specimens were correctly identified 100% of the time.

BIBLIOGRAPHY

- 1. Clémence Magnard, Martine Valette, Michèle Aymard, Bruno Lina. Comparison of Two Nested PCR, Cell Culture, and Antigen Detection for the Diagnosis of Upper Respiratory Tract Infections due to Influenza Viruses. Journal of Medical Virology, 59:215–220 (1999)
- 2. Shuo Su, Gary Wong, Weifeng Shi, et al. Epidemiology, Genetic recombination, and pathogenesis of coronaviruses. Trends in Microbiology, June 2016, vol. 24, No. 6: 490-502
- 3. Susan R. Weiss, Julian L. Leibowitz, Coronavirus Pathogenesis, Advances in Virus Research, Volume 81: 85-164

Index of Symbols	
	Manufacturer
	<i>In vitro</i> diagnostic medical device
	Consult instructions for use
	Authorized representative in the European Community
	Contains sufficient for <n> tests
	Use-by date
	Batch code
	Medical device
	Temperature limit
	Do not re-use
	Catalogue number
	Date of manufacture

Index of Contents

Negative Control Swab	Negative Control Swab
Positive Control Swab	Positive Control Swab
Extraction Buffer Tubes	Extraction Buffer Tubes
Disposable Swabs	Disposable Swabs
Nasal Swabs	Nasal Swabs
Nasopharyngeal Swabs	Nasopharyngeal Swabs

**ACON Biotech (Hangzhou) Co., Ltd.**
No.210 Zhenzhong Road, West Lake District, Hangzhou, P.R.China, 310030

**MedNet EC-REP GmbH**
Borkstrasse 10
48163 Muenster, Germany

Number: 1151485401
Effective date: 2022-09-06