JOY TEST

Lyme IgG/IgM Rapid Test Cassette (Whole Blood/Serum/Plasma) Lymská borelióza test (IgG/IgM) (Plná krev/sérum/plazma)

Návod

REF ILY-402

Čeština

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Lyme IgG/IgM Rapid Test Cassette (Whole Blood/Serum/Plasma) Package Insert

REF ILY-402
English

A rapid test for the qualitative detection of IgG and IgM antibodies to *Borrelia* in human whole blood, serum or plasma specimens.

For professional *in vitro* diagnostic use only.

INTENDED USE

The Lyme IgG/IgM Rapid Test Cassette is a lateral flow chromatographic immunoassay for the qualitative detection of IgG and IgM antibodies to *Borrelia* in human whole blood, serum or plasma specimen.

SUMMARY

Lyme disease, also known as Lyme borreliosis, is an infectious disease caused by bacteria of the *Borrelia* sp. which is spread by ticks.² The most common sign of infection is an expanding area of redness on the skin, known as erythema migrans, that begins at the site of a tick bite about a week after it has occurred.¹ The rash is typically neither itchy nor painful. Approximately 25–50% of infected people do not develop a rash.¹ Other early symptoms may include fever, headache and feeling tired.¹ If untreated, symptoms may include loss of the ability to move one or both sides of the face, joint pains, severe headaches with neck stiffness, or heart palpitations, among others.¹ Months to years later, repeated episodes of joint pain and swelling may occur.¹ Occasionally, people develop shooting pains or tingling in their arms and legs.¹ Despite appropriate treatment, about 10 to 20% of people develop joint pains, memory problems, and feel tired for at least six months.^{1,4}

Lyme disease is transmitted to humans by the bite of infected ticks of the genus *Ixodes*.⁵ Usually, the tick must be attached for 36 to 48 hours before the bacteria can spread.⁶ In North America, *Borrelia burgdorferi* and *Borrelia mayonii* are the causes.^{2,7} In Europe and Asia, the bacteria *Borrelia afzelii* and *Borrelia garinii* are also causes of the disease.² The disease does not appear to be transmissible between people, by other animals, or through food.⁶ Diagnosis is based upon a combination of symptoms, history of tick exposure, and possibly testing for specific antibodies in the blood.^{3,8} Blood tests are often negative in the early stages of the disease.² Testing of individual ticks is not typically useful.⁹

PRINCIPLE

The Lyme IgG/IgM Rapid Test Cassette (Whole Blood/Serum/Plasma) is a qualitative membrane-based immunoassay for the detection of IgG and IgM antibodies to *Borrelia* in whole blood, serum or plasma specimens. This test consists of two components, an IgG component and an IgM component. In the IgG component, anti-human IgG is coated in IgG test line region. During testing, the specimen reacts with *Borrelia* antigen-coated particles in the test cassette. The mixture then migrates upward on the membrane chromatographically by capillary action and reacts with the anti-human IgG in IgG test line region, if the specimen contains IgG antibodies to *Borrelia*. A colored line will appear in IgG test line region as a result of this. Similarly, anti-human IgM is coated in IgM test line region and if specimen contains IgM antibodies to *Borrelia*, the conjugate-specimen complex reacts with anti-human IgM. A colored line will appear in IgM test line region as a result.

Therefore, if the specimen contains anti-*Borrelia* IgG antibodies, a colored line will appear in IgG test line region. If the specimen contains anti-*Borrelia* IgM antibodies, a colored line will appear in IgM test line region. If the specimen does not contain anti-*Borrelia* antibodies, no colored line will appear in either of the test line regions, indicating a negative result. To serve as a procedural control, a colored line will always appear in the control line region, indicating that the proper volume of specimen has been added and membrane wicking has occurred.

REAGENTS

The test contains anti-human IgM and anti-human IgG as the capture reagent, *Borrelia* antigen as the detection reagent. A goat anti-human IgG is employed in the control line system.

PRECAUTIONS

- For professional *in vitro* diagnostic use only. Do not use after expiration date.
- Do not eat, drink or smoke in the area where the specimens or kits are handled.
- Do not use test if pouch is damaged.
- Handle all specimens as if they contain infectious agents. Observe established precautions against microbiological hazards throughout all procedures and follow the standard procedures for proper disposal of specimens.
- Wear protective clothing such as laboratory coats, disposable gloves and eye protection when specimens are assayed.
- The used test should be discarded according to local regulations.
- Humidity and temperature may adversely affect results.

STORAGE AND STABILITY

Store as packaged in the sealed pouch at room temperature or refrigerated (2-30°C). The test is stable through the expiration date printed on the sealed pouch. The test must remain in the sealed pouch until use. **DO NOT FREEZE.** Do not use beyond the expiration date.

SPECIMEN COLLECTION AND PREPARATION

- The Lyme IgG/IgM Rapid Test Cassette (Whole Blood/Serum/Plasma) can be performed using whole blood (from venipuncture or fingerstick), serum or plasma.
- To collect **Fingerstick Whole Blood Specimens**:
 - Wash the patient's hand with soap and warm water or clean with an alcohol swab. Allow to dry.
 - Massage the hand without touching the puncture site by rubbing down the hand towards the fingertip of the middle or ring finger.
 - Puncture the skin with a sterile lancet. Wipe away the first sign of blood.
 - Gently rub the hand from wrist to palm to finger to form a rounded drop of blood over the puncture site.
- Add the Fingerstick Whole Blood specimen to the test by using **a capillary tube**:
 - Touch the end of the capillary tube to the blood until filled to approximately 10µL. Avoid air bubbles.
 - Place the bulb onto the top end of the capillary tube, then squeeze the bulb to dispense the whole blood to the specimen area of the test cassette.
- Separate serum or plasma from blood as soon as possible to avoid hemolysis. Use only clear non-hemolyzed specimens.
- Testing should be performed immediately after the specimens have been collected. Do not leave the specimens at room temperature for prolonged periods. Serum and plasma specimens may be stored at 2-8°C for up to 3 days, for long term storage, specimens should be kept below -20°C. Whole blood collected by venipuncture should be stored at 2-8°C if the test is to be run within 2 days of collection. Do not freeze whole blood specimens. Whole blood collected by fingerstick should be tested immediately.
- Bring specimens to room temperature prior to testing. Frozen specimens must be completely thawed and mixed well prior to testing. Specimens should not be frozen and thawed repeatedly.
- If specimens are to be shipped, they should be packed in compliance with local regulations covering the transportation of etiological agents.
- EDTA K2, Heparin sodium, Citrate sodium and Potassium Oxalate can be used as the anticoagulant for collecting the specimen.

MATERIALS

Materials provided

- Test cassettes
- Droppers
- Package insert
- Buffer

Materials required but not provided


- Specimen collection containers
- Centrifuge
- Lancets (for fingerstick whole blood only)
- Timer
- Pipette and disposable tips (optional)

DIRECTIONS FOR USE

Allow the test, specimen, buffer and/or controls to reach room temperature (15-30°C) prior to testing.

- Bring the pouch to room temperature before opening it. Remove the test cassette from the sealed pouch and use it as soon as possible.
- Place the cassette on a clean and level surface.
For **Serum/Plasma** specimen:
 - Use a **pipette**: To transfer **5 µL of Serum/Plasma** to the specimen well (S), then add **3 drops of buffer (approximately 120 µL) to the buffer well (B)**.
 - Use a **dropper**: Hold the dropper vertically, draw the specimen up to the **upper end of the nozzle** as shown in illustration below (**approximately 5 µL**). Transfer the specimen to the specimen well (S), then add **3 drops of buffer (approximately 120 µL) to the buffer well (B)**, and start the timer.
- For **Whole Blood** specimen:
 - Use a **pipette**: To transfer **10 µL of whole blood** to the specimen well (S), then add **3 drops of buffer (approximately 120 µL) to the buffer well (B)**.
 - Use a **dropper**: Hold the dropper vertically, draw the specimen about **1 cm above the upper end of the nozzle** and transfer **1 full drop (approx. 10 µL)** of specimen to the sample well (S). Then add **3 drops of buffer (approximately 120 µL) to the buffer well (B)**, and start the timer.
- Wait for the colored line(s) to appear. Read results at **10 minutes**. Do not interpret the result after **20 minutes**.

Note: It is suggested not to use the buffer beyond 3 months after opening the vial.

INTERPRETATION OF RESULTS

IgG POSITIVE: * **Two colored lines appear.** One colored line should always appear in the control line region (C) and another line should be in the IgG line region.

IgM POSITIVE: * **Two colored lines appear.** One colored line should always appear in the control line region (C) and another line should be in the IgM line region.

IgG and IgM POSITIVE: * **Three colored lines appear.** One colored line should always appear in the control line region (C) and two test lines should be in the IgG line region and IgM line region.

***NOTE:** The intensity of the color in the test line regions may vary depending on the concentration of anti-Lyme antibodies present in the specimen. Therefore, any shade of color in the test line region should be considered positive.

NEGATIVE: **One colored line appears in the control line region (C).** No line appears in the IgG region or the IgM region.

INVALID: Control line fails to appear. Insufficient specimen volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test with a new test. If the problem persists, discontinue using the test kit immediately and contact your local distributor.

QUALITY CONTROL

Internal procedural controls are included in the test. A colored line appearing in the control region (C) is an internal procedural control. It confirms sufficient specimen volume and correct procedural technique. Control standards are not supplied with this kit; however, it is recommended that positive and negative controls be tested as good laboratory testing practice to confirm the test procedure and to verify proper test performance.

LIMITATIONS

- The Lyme IgG/IgM Rapid Test Cassette (Whole Blood/Serum/Plasma) is for *in vitro* diagnostic use only. This test should be used for detection of IgG and IgM antibodies to *Borrelia* in whole blood, serum or plasma specimens. Neither the quantitative value nor the rate of increase in the concentration of IgG or IgM antibodies to *Borrelia* can be determined by this qualitative test.
- The Lyme IgG/IgM Rapid Test Cassette (Whole blood/Serum/Plasma) will only indicate the presence of IgG and IgM antibodies to *Borrelia* in the specimen and should not be used as the sole criteria for the diagnosis of Lyme infections.
- As with all diagnostic tests, all results must be considered with other clinical information available to the physician.
- If the test result is negative and clinical symptoms persist, additional follow-up testing using other clinical methods is suggested. A negative result at any time does not preclude the possibility of *Borrelia* infection.
- The hematocrit level of the whole blood can affect the test results.
- Hematocrit level needs to be between 25% and 65% for accurate results.

PERFORMANCE CHARACTERISTICS

Sensitivity and Specificity

The Lyme IgG/IgM Rapid Test Cassette (Whole Blood/Serum/Plasma) was compared with a leading commercial ELISA Lyme IgG tests and ELISA Lyme IgM tests; the results show that Lyme IgG/IgM Rapid Test Cassette (Whole Blood/Serum/Plasma) has a high sensitivity and specificity.

IgG Results

Lyme IgG/IgM Rapid Test for IgG	Method	ELISA		Total Results
	Results	Positive	Negative	
	Positive	21	1	22
Total Results		1	89	90
		22	90	112

Relative Sensitivity: 95.5% (95%CI*: 87.3%-100%)

Relative Specificity: 98.9% (95%CI*: 97.1%-99.8%)

Accuracy: 98.2% (95%CI*: 93.7%-99.8%)

*Confidence Interval

IgM Results

Lyme IgG/IgM Rapid Test for IgM	Method	ELISA		Total Results
	Results	Positive	Negative	
	Positive	17	1	18
Total Results		1	89	90
		18	90	108

Relative Sensitivity: 94.4% (95%CI*: 72.7%-99.9%)

Interval

Relative Specificity: 98.9% (95%CI*: 96.7%-100%)

Accuracy: 98.1% (95%CI*: 93.5%-99.8%)

*Confidence

Precision Intra-Assay

Within-run precision has been determined by using 3 replicates of five specimens: negative, IgG low positive, IgG high positive, IgM low positive, IgM high positive. The negative, low positive, and high positive values were correctly identified >99% of the time.

Inter-Assay

Between-run precision has been determined by 3 independent assays on the same specimens: negative, IgG low positive, IgG high positive, IgM low positive, IgM high. Three different lots of the Lyme IgG/IgM Rapid Test cassette (Whole Blood/Serum/Plasma) have been tested over a 3-days period using negative, low positive, and high positive specimens. The specimens were correctly identified >99% of the time.

Cross-reactivity

The Lyme IgG/IgM Rapid Test Cassette (Whole Blood/Serum/Plasma) has been tested for anti-HAV IgM, HBsAg, anti-HCV IgG, anti-HIV IgG, anti-RF IgG, anti-Syphilis IgG, anti-*H. Pylori* IgG, anti-Rubella IgG, anti-Toxo IgG, anti-HSV 1 IgG, anti-HSV 2 IgG , anti-CMV IgG, anti-Rubella IgM, anti-Toxo IgM, anti-HSV 1 IgM, anti-HSV 2 IgM and anti-CMV IgM positive specimens. The results showed no cross-reactivity.

Interfering Substances

The following compounds have been tested using the LymeIgG/IgM Rapid Test Cassette (Whole Blood/Serum/Plasma) and no interference was observed.

Acetaminophen: 20 mg/dL	Caffeine: 20 mg/dL	Acetylsalicylic Acid: 20 mg/dL
Gentisic Acid: 20 mg/dL	Ascorbic Acid: 2 g/dL	Albumin: 2 g/dL
Creatin: 200 mg/dL	Hemoglobin 1000 mg/dL	Bilirubin: 1 g/dL
Oxalic Acid: 60 mg/dL		

BIBLIOGRAPHY

1. "Signs and Symptoms of Lyme Disease".cdc.gov. 11 January 2013. Archived from the original on 16 January 2013. Retrieved 2 March 2015.
2. Shapiro, ED (1 May 2014). "Clinical practice. Lyme disease" (PDF). N. Engl. J. Med. 370 (18): 1724–31.
3. "Lyme Disease Diagnosis and Testing". cdc. gov. 10 January 2013. Archived from the original on 2 March 2015. Retrieved 2 March 2015.
4. Aucott JN (2015). "Posttreatment Lyme disease syndrome". Infect. Dis. Clin. N. Am. 29 (2): 309–23.
5. Johnson RC (1996). "Borrelia". In Baron S; et al. Baron's Medical Microbiology (4th ed.). Univ of Texas Medical Branch. ISBN0-9631172-1-1. PMID21413339. Archived from the original on 7 February 2009.
6. "Lyme disease transmission". cdc. gov. 11 January 2013. Archived from the original on 3 March 2015. Retrieved 2 March 2015.
7. Jump up^Pritt, BS; Mead, PS; Johnson, DK; Neitzel, DF; RespicioKingry, LB; Davis, JP; Schiffman, E; Sloan, LM; Schriefer, ME; Replogle, AJ; Paskewitz, SM; Ray, JA; Bjork, J; Steward, CR; Deedon, A; Lee, X; Kingry, LC; Miller, TK; Feist, MA; Theel, ES; Patel, R; Irish, CL; Petersen, JM (5 February 2016). "Identification of a novel pathogenic Borrelia species causing Lyme borreliosis with unusually high spirochaetaemia: a descriptive study". Lancet Infect. Dis. 16: 556–564.
8. "Two-step Laboratory Testing Process". cdc. gov. 15 November 2011. Archived from the original on 12 March 2015. Retrieved 2 March 2015.
9. Jump up^"Testing of Ticks".cdc.gov. 4 June 2013.Archivedfrom the original on 19 February 2015. Retrieved2 March2015.

	Caution
	For <i>in vitro</i> diagnostic use only
 2°C / 30°C	Store between 2-30°C
	Do not use if package is damaged
	Authorized representative in EU
	Catalog #
	Tests per kit
	Use by
	Lot number
	Manufacturer
	Do not reuse
	Consult instructions for use
	Importer
	Distributor

Manufacturer

Hangzhou Alltest Biotech Co., Ltd.
#550, Yinhai Street
Hangzhou Economic & Technological Development Area
Hangzhou, 310018 P.R. China
www.alltests.com.cn Email: info@alltests.com.cn

MedNet EC-REP GmbH
Borkstrasse 10
48163 Muenster
Germany
Czech Original Products s.r.o.
Koulova 6, Praha 6, 160 00 - CZ
JOYMED.cz - ID: 08595771
Czech Original Products s.r.o.
Koulova 6, Praha 6, 160 00 - CZ
JOYMED.cz - ID: 08595771

Number:
Revision date: 2025-06-05