JOY

TEST

GISTA-702H(1T/Kit);GISTA-702H(2T/Kit);GISTA-702H(5T/Kit);GISTA-702H(25T/Kit)

English

Rapid test for the qualitative detection of group A *Streptococcus* antigens in human throat swab samples. For self-testing, *in-vitro* diagnostic use only.

INTENDED PURPOSE

The Strep A Rapid Test is a lateral flow chromatographic immunoassay for the qualitative detection of group A *Streptococcus* antigens in human throat swab samples. The test is not automated. This product is suitable for self-testing by patients over 16 years of age. It is recommended that individuals aged from 3-16 should be tested by a parent or legal guardian. The test is an *in-vitro* diagnostic device intended for use as an aid in the diagnosis of a *Streptococcus* A (Strep A) infection in patients with typical symptoms, such as fever, pharyngitis, etc.

SUMMARY

Streptococcus pyogenes is a species of non - motile gram - positive cocci which are categorised as Lancefield group A, which may cause serious infections such as pharyngitis, respiratory infection, impetigo, endocarditis, meningitis, puerperal sepsis and arthritis.¹ Left untreated, these infections may lead to serious complications, including rheumatic fever, peritonsillar abscesses or scarlet fever with a rash and sore throat.² Traditional identification procedures for group A streptococcal infection involve the isolation and identification of viable organisms using techniques that require 24 to 48 hours or longer.^{3,4}

The Strep A Rapid Test is an immunological test designed to specifically identify *Streptococcus* A bacteria in oropharyngeal infections. The test helps the patient understand whether a sore throat is caused by *Streptococcus* A.

TEST PRINCIPLE

The Strep A Rapid Test is a lateral flow immunoassay based on the principle of the double-antibody sandwich technique. The test cassette consists of: 1) a burgundy conjugate pad containing rabbit anti-group A *Streptococcus* monoclonal antibodies conjugated with colloidal gold (Strep A Ab conjugates), 2) a nitrocellulose membrane strip containing a test line (T line) and a control line (C line). The T line is pre-coated with further rabbit anti-group A *Streptococcus* monoclonal antibodies, and the C line is pre-coated with goat anti-rabbit IgG.

When an adequate volume of a test specimen is dispensed into the sample well of the test cassette, the specimen migrates by capillary action along the membrane. The group A *Streptococcus* antigens, if present in the specimen, will bind to the Strep A Ab conjugates. The immunocomplex is then captured on the membrane by the pre-coated rabbit anti-group A *Streptococcus* monoclonal antibody reagent, forming a burgundy T line, indicating a group A *Streptococcus* positive test result.

The test contains an internal control (C line) which, regardless of colour development on the T line, should exhibit a burgundy coloured control line, caused by the immunocomplex of the control antibodies . If no C line appears, the test result is invalid and the specimen must be retested using another test cassette.

WARNINGS AND PRECAUTIONS

- For *in-vitro* diagnostic use only.
- Read these instructions for use carefully before performing the test. The test cassette is reliable if the instructions are followed correctly (reaction time, how the sample is collected, sample extraction etc.).
- The assay is intended for use as an aid in the diagnosis of a group A streptococcal infection. The test does not determine the state of infection or associated diseases.
- Dispose of following use. The test cassette is not suitable for reuse.
- Keep the test out of the reach of children.
- Extraction Reagent 1 contains sodium nitrite. Warning: harmful if swallowed, causes serious eye irritation.
- Do not touch the reaction area on the test cassette.
- Do not use beyond the expiry date.
- Do not use the test cassette if the pouch is punctured or poorly sealed.
- The test is for external use only. Do not swallow reagents 1 and 2 or other items in the box. If swallowed, seek medical advice immediately and show the remaining contents of the box, the instructions for use and the packaging.
- If liquids come into contact with the eyes, skin or other mucosae, immediately rinse thoroughly with water and contact a physician. Show the bottles' labels.
- Do not interchange bottle caps between reagents.
- Please follow local regulations when disposing of used tests.
- All test kit components are for single-use only, with the exception of the tube stand (which can be re-used), the instructions for use and reagents 1 and 2, which must be stored tightly sealed, at a temperature between 2° C and 30° C. The rest of the test kit should be stored at 2 ~ 30° C.
- The swab is a single-use sterilised medical device and is provided in a sterilised state. It has been sterilised using ethylene oxide, and CANNOT be re-used. Do not use if the packaging is damaged.
- If any serious events occur in relation to the test, please report them to Greylynx and the competent authority of your country.
- In the event that the patient is unable to complete the test themselves, it is recommended that the test be performed with the help of an adult. If you are currently helping someone else with the test, please follow the instructions for completing the test.

COMPONENTS

The Strep A Rapid Test contains the 'test cassette (packaged in pouch with desiccant)', 'disposable swab', 'tongue depressor', 'Reagent 1', 'Reagent 2', 'test tube', 'tube tip', 'tube stand' and 'Instructions for use'.

Test cassette

Disposable swab

Tongue depressor

Reagent 1&2

Test tube & tip

Tube stand

Instructions for use

• **Test cassette:** The test cassette is composed of glass fiber, a nitrocellulose membrane, a plastic backer, absorbent paper. The main components are rabbit anti-group A *Streptococcus* monoclonal antibodies, goat anti-rabbit polyclonal antibodies, colloidal gold labeled rabbit anti-group A *Streptococcus* monoclonal antibodies.

• **Disposable swab:** For sample collection.

• **Reagent 1:** Bottle with 1 M sodium nitrite, with white bottle cap

• **Reagent 2:** Bottle with 0.15 M acetic acid, with red bottle cap

• **Test tube:** For mixing reagent 1, reagent 2 and processing samples

• **Tube tip:** To cover the test tube, drip the treated sample

• **Tube stand:** Holds the test tube (Note: Please use the opening marked on the test kit box as a tube stand.)

• **Tongue depressor:** Fixate the tongue to assist in specimen sampling

• **Instructions for use**

Components of the Strep A Rapid Test

Package

GISTA-702H(1T/Kit)

GISTA-702H(2T/Kit)

GISTA-702H(5T/Kit)

GISTA-702H(25T/Kit)

1 test/kit

2 tests/kit

5 tests/kit

25 tests/kit

Test cassette

Disposable swab

Reagent 1

Reagent 2

Test tube

Tube tip

Tube stand

Tongue depressor

Instructions for use

MATERIALS REQUIRED BUT NOT SUPPLIED

• Timer

STORAGE AND STABILITY

- Store the test kit between 2-30°C and 10% to 80 % relative humidity. Keep away from light. Exposure to temperatures and / or humidity outside the specified conditions may cause inaccurate results.
- Do not freeze. Use the test at temperatures between 15-30°C.
- Do not use beyond the expiration date (printed on the foil pouch and box).
- To avoid any risk and obtain optimal results, the test should be performed within one hour of removing the test cassette from its sealed foil pouch. When the temperature is higher than 30°C or the humidity is higher than 90%, use immediately after opening the foil pouch.
- The reagents are stable for 24 months. Extraction reagents, stressed by multiple cycles of opening/closing, still demonstrate the expected functionality.

Note: All expiration dates are printed in Year-Month-Day format. 2022-06-18 indicates June 18, 2022.

SPECIMEN COLLECTION AND PREPARATION

- Specimen collection**
 - Collect a throat swab sample using the disposable swab provided with the test kit. Swab the posterior pharynx, tonsils, and other surrounding inflamed areas, avoiding contact with the tongue, inner cheeks, and teeth. Fasting more than half an hour before collecting a sample is recommended to avoid food interference.
- Specimen handling**
 - Testing is recommended immediately after collecting a throat swab sample. If testing cannot be performed immediately, throat swab samples should be stored in a clean, dry, airtight container at room temperature for up to 4 hours and at 2-8°C for up to 24 hours. DO NOT FREEZE.

TEST PROCEDURE

Allow the tests, samples, and reagents to equilibrate to room temperature(15-30°C) prior to testing.

- If the specific patient age is 3-16, the following test steps should be carried out by a parent or legal guardian.
- Wash hands with hot water and soap, rinse with clean water and dry.
- Place the tube stand on a clean, dry, flat surface. Place the empty plastic test tube in one of the openings in the tube stand-FIG.1

FIG-1

4. Add 6 drops of Reagent 1 (white bottle cap) and 4 drops of Reagent 2 (red bottle cap) to the test tube, and fully mix.-FIG.2, FIG.3

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