



Ref.	Size	Unit	Box Dimensions	EAN Box	Carton Dimensions	EAN Carton	Qty per pallet
655130	XS	10 x 100	230*125*60mm	8720598882231	315*258*245mm	8720598882148	80 cartons / 800 boxes
655131	S	10 x 100	230*125*60mm	8720598882255	315*258*245mm	8720598882162	80 cartons / 800 boxes
655132	M	10 x 100	230*125*60mm	8720598882279	315*258*245mm	8720598882186	80 cartons / 800 boxes
655133	L	10 x 100	230*125*60mm	8720598882293	315*258*245mm	8720598882209	80 cartons / 800 boxes
655134	XL	10 x 100	230*125*60mm	8720598882316	315*258*245mm	8720598882223	80 cartons / 800 boxes



STERISETS INTERNATIONAL B.V.  
Ketelmeer 3 | 5347 JX Oss | The Netherlands  
[www.sterisets.com](http://www.sterisets.com)

## Features:

- Latex, powder and odor free.
- Excellent softness and wearing fitness.
- Braided cuff ensures easy donning and prevents rolling down.
- Finger textured for a reliable, consistent grip.
- 600% stretch.
- 5 year' shelf life

## Quality Standards:

- EN455-1:2020, EN455-2:2015, EN455-3:2015, EN455-4:2009
- EN ISO 374-1:2016+A1:2018
- EN ISO 374-5:2016
- EN ISO 21420:2020
- AQL <1.5

## Applications:

- Medical Purpose / Examination
- Healthcare and nursing
- Industrial purpose / PPE
- General housekeeping
- Laboratory
- IT Industry